




Pockels effect


- An anisotropic crystal such as CZT will have different refractive indices based on crystal lattice orientation. This phenomenon is also called “birefringence”.
- Pockels effect is an electro-optical effect of *varying the directionally dependent refractive index* in response to an applied E-field.
 - E-Field “modulates” the refractive indices in the crystal
- Most application use Pockels effect to change the state of light.
 - Optical shutter, attenuator, ultrafast lasers
- We use the resultant light to investigate the E-field profile of the material



Animation of Pockels cell



Pockels apparatus



- First polarizer is fixed at 45 degrees
- Second polarizer (aka “analyzer”) has two orientations:
 - Parallel (45 degrees) for calibration
 - Crossed (135 degrees) for bias measurements


$$T = \frac{I_{\perp,bias} - I_{\perp,0V}}{I_{\parallel,0V} - I_{\parallel,bkg}} = \sin^2(\alpha E)$$

- $$E = \frac{1}{\alpha} \sqrt{\arcsin(T)}$$

$$\left(\alpha = \frac{\sqrt{3} \pi n_0^3 r d}{2 \lambda_0} \right)$$
- Both $T(x, y)$ and $E(x, y)$ are 2D maps
- When there’s no electric bias (0V), $T \approx 0$


Calculation of E-field map

- $T(x, y) = \frac{I_{\perp, bias} - I_{\perp, 0V}}{I_{\parallel, 0V} - I_{\parallel, bkg}} = \sin^2(\alpha E(x, y))$
- Capture 3 calibration images
 - Light on, analyzer at parallel, $I_{\parallel, 0V}$
 - Light on, analyzer at crossed, $I_{\perp, 0V}$
 - Light off, analyzer at parallel, $I_{\parallel, bkg}$
- Pockels effect is induced for any electric bias, start to see transmission of light
- E-field map is calculated by inverting the equation above



E-field profile per applied voltage

- The total E-field within CZT increases with bias
- The slope of each curve tells us the space charge density inside



How to calculate space charge density

Gauss's Law

$$\oint_S \vec{E} \cdot d\vec{A} = \frac{Q_{enclosed}}{\epsilon_0} = \frac{\rho \cdot V}{\epsilon_0}$$

Assume an infinite slab of uniform charge density, the internal electric field is

$$E(x) = \frac{\rho}{\epsilon_0} \left(\frac{d}{2} - x \right)$$

$$E_{total}(x) = E(x) + E_{bias}$$

